

Doing quantization and getting overall metric for all the models

```

import tensorflow as tf
import pydicom
import numpy as np
import matplotlib.pyplot as plt
import os
import pandas as pd
import glob
from tqdm import tqdm
import seaborn as sns
import warnings
warnings.filterwarnings("ignore")
import gc
gc.collect()
import cv2
import tensorflow as tf
import keras
from PIL import Image, ImageDraw
from PIL import ImagePath
import imgaug.augmenters as iaa
from skimage import exposure
from tensorflow import keras
from tensorflow.keras.layers import *
from tensorflow.keras.preprocessing import image
from tensorflow.keras.models import Model, load_model
from tensorflow.keras.layers import UpSampling2D
from tensorflow.keras.layers import MaxPooling2D, GlobalAveragePooling2D
from tensorflow.keras.layers import concatenate
from tensorflow.keras.layers import Multiply
from tensorflow.keras.callbacks import EarlyStopping, ModelCheckpoint
from tensorflow.keras import backend as K
from tensorflow.keras.layers import Input, Add, Dense, Activation, ZeroPadding2D, BatchNormal
from tensorflow.keras.models import Model, load_model
from tensorflow.keras.utils import plot_model
from tensorflow.keras.initializers import glorot_uniform, he_normal
from tensorflow.keras.losses import binary_crossentropy
import pathlib
K.set_image_data_format('channels_last')

import numpy as np

def mask2rle(img, width, height):
    rle = []
    lastColor = 0;
    currentPixel = 0;
    runStart = -1;
    runLength = 0;

```

```

runLength = 0,

for x in range(width):
    for y in range(height):
        currentColor = img[x][y]
        if currentColor != lastColor:
            if currentColor == 255:
                runStart = currentPixel;
                runLength = 1;
            else:
                rle.append(str(runStart));
                rle.append(str(runLength));
                runStart = -1;
                runLength = 0;
                currentPixel = 0;
        elif runStart > -1:
            runLength += 1
        lastColor = currentColor;
        currentPixel+=1;

return " ".join(rle)

def rle2mask(rle, width, height):
    mask= np.zeros(width* height)
    array = np.asarray([int(x) for x in rle.split()])
    starts = array[0::2]
    lengths = array[1::2]

    current_position = 0
    for index, start in enumerate(starts):
        current_position += start
        mask[current_position:current_position+lengths[index]] = 255
        current_position += lengths[index]

    return mask.reshape(width, height)

#https://www.kaggle.com/jesperdramsch/intro-chest-xray-dicom-viz-u-nets-full-data
def getInfoDICOM(path, data_rle_train, toPrint = False, train = True):

    info = {}

    path = os.path.join(path)
    #reading the data using methods of pydicom
    data = pydicom.dcmread(path)
    info['path'] = path
    info['age'] = data.PatientAge
    info['sex']= data.PatientSex
    info['ImageId'] = data.SOPInstanceUID

    if train: #test doesn't have encoded pixels data, we have to predict them
        #SOPInstanceUID contains the storage type which resembles the image id given in train

```

```
#SOPInstanceUID contains the storage type which resembles the image ids given in train
encodedPixels = data_rle_train[data_rle_train['ImageId'] == data.SOPInstanceUID]["EncodedPixels"]
info['encodedPixels'] = encodedPixels
```

```
info['lenOfEncodedPixels'] = len(encodedPixels)
```

```
#this is for visualization purpose
```

```
if '-1' in encodedPixels or len(encodedPixels) == 0:
```

```
    info['has_pneumothorax'] = 0
```

```
else:
```

```
    info['has_pneumothorax'] = 1
```

```
if toPrint:
```

```
    print("Path.....:", path)
```

```
    print("Patient's Name.....:", data.PatientName)
```

```
    print("Patient's Id.....:", data.PatientID)
```

```
    print("Patient's Age.....:", data.PatientAge)
```

```
    print("Patient's Sex.....:", data.PatientSex)
```

```
    print("Original X Ray")
```

```
    plt.figure(figsize=(10,10))
```

```
    plt.imshow(data.pixel_array, cmap=plt.cm.bone)
```

```
    plt.show()
```

```
return info
```

```
def computeMasks(train, path):
```

```
    mask_paths = []
```

```
    for index, row in train.iterrows():
```

```
        pixels = row['encodedPixels']
```

```
        mask = np.zeros((1024, 1024))
```

```
        if row['has_pneumothorax'] == 1:
```

```
            for pix in pixels:
```

```
                mask = mask + rle2mask(pix, 1024, 1024).T
```

```
        img = Image.fromarray(mask).convert('L')
```

```
        path2save = path + str(index+1) + ".jpg"
```

```
        img.save(path2save)
```

```
        mask_paths.append(path2save)
```

```
    train["mask"] = mask_paths
```

```
    return train
```

```
def convertImagesToJpeg(train, path):
```

```
    images_paths = []
```

```
    for index, row in train.iterrows():
```

```
        path2dicom = row['path']
```

```
        ds = pydicom.read_file(path2dicom) # read dicom image
```

```
        img = ds.pixel_array # get image array
```

```
        img_mem = Image.fromarray(img) # Creates an image memory from an object exporting the
```

```
        path2save = path + str(index + 1) + ".jpg"
```

```
        img_mem.save(path2save)
```

```

img_name.save(path2save,
              images_paths.append(path2save)
train["images_paths"] = images_paths
return train

```

```
img_size = 256
```

```

def load_data():
    print("Loading Data ...")
    data_rle_train = pd.read_csv("./siim/train-rle.csv")
    #The second column contains a space infront of it so manually rewriting it
    data_rle_train.columns = ['ImageId', 'EncodedPixels']

    print("Extracting info from dicom")
    filesList = glob.glob('./siim/dicom-images-train/**/*.dcm')
    train = pd.DataFrame()
    train_list = []
    for file in filesList:
        info = getInfoDICOM(file, data_rle_train, False, True)
        if info['has_pneumothorax'] == 1:
            train_list.append(info)
    train = pd.DataFrame(train_list)
    train.loc[train.lenOfEncodedPixels == 0 , 'encodedPixels'] = np.array(-1)
    print("Extracted info from dicom")

    print("Computing Masks")
    path = "siim/dicom-mask-train/"
    try:
        os.makedirs(path)
    except:
        pass
    train = computeMasks(train, path)
    print("Computed Masks")

    print("Converting Images to Jpg")
    path = "siim/dicom-images-train-jpg/"
    try:
        os.makedirs(path)
    except:
        pass
    train = convertImagesToJpeg(train, path)
    print("Converted Images to Jpg")

    return train

threshold = 0.38 #learned from previous model
train = load_data()

```

```

Loading Data ...
Extracting info from dicom

```

Extracted info from dicom
 Computing Masks
 Computed Masks
 Converting Images to Jpg
 Converted Images to Jpg

```
train.head()
```

| | path | age | sex | | ImageId | enco |
|---|---|-----|-----|---|---------|--------------|
| 0 | ./siim/dicom-images-train/1.2.276.0.7230010.3.... | 26 | M | 1.2.276.0.7230010.3.1.4.8323329.1051.151787516... | | [120: 37 62] |
| 1 | ./siim/dicom-images-train/1.2.276.0.7230010.3.... | 16 | M | 1.2.276.0.7230010.3.1.4.8323329.4168.151787518... | | [166: 1 9] |
| 2 | ./siim/dicom-images-train/1.2.276.0.7230010.3.... | 63 | F | 1.2.276.0.7230010.3.1.4.8323329.5057.151787518... | | [672: 8 10] |
| 3 | ./siim/dicom-images-train/1.2.276.0.7230010.3.... | 13 | F | 1.2.276.0.7230010.3.1.4.8323329.3444.151787517... | | 10 23 9 |
| 4 | ./siim/dicom-images- | 51 | F | 1.2.276.0.7230010.3.1.4.8323329.1478.151787516 | | [26: 45] |

```
smooth = 1e-5
```

```
def dice_coef(y_true, y_pred):
    y_true_f = tf.keras.layers.Flatten()(y_true)
    y_pred_f = tf.keras.layers.Flatten()(y_pred)
    intersection = tf.reduce_sum(y_true_f * y_pred_f)
    return (2. * intersection + smooth) / (tf.reduce_sum(y_true_f) + tf.reduce_sum(y_pred_f))
```

```
def combined_loss(y_true, y_pred):
    return 1.0 - dice_coef(y_true, y_pred) + binary_crossentropy(y_true, y_pred)
```

```
!unzip pneumothorax_predictor.zip
```

```
Archive: pneumothorax_predictor.zip
  creating: pneumothorax_predictor/pneumothorax_predictor/
  creating: pneumothorax_predictor/pneumothorax_predictor/assets/
 inflating: pneumothorax_predictor/pneumothorax_predictor/saved_model.pb
  creating: pneumothorax_predictor/pneumothorax_predictor/variables/
 inflating: pneumothorax_predictor/pneumothorax_predictor/variables/variables.data-0000
 inflating: pneumothorax_predictor/pneumothorax_predictor/variables/variables.index
```

```
pneumothorax_predictor = tf.keras.models.load_model("pneumothorax_predictor/pneumothorax_pred
```

```
#Predictor
```

```
converter_pred = tf.lite.TFLiteConverter.from_keras_model(pneumothorax_predictor)
```

```
tflite_models_dir = pathlib.Path("/content/pneumothorax/")
```

```
tflite_models_dir.mkdir(exist_ok=True, parents=True)
```

```
converter_pred.optimizations = [tf.lite.Optimize.DEFAULT]
```

```
tflite_predict_model = converter_pred.convert()
```

```
tflite_model_quant_file = tflite_models_dir/"predictor.tflite"
```

```
tflite_model_quant_file.write_bytes(tflite_predict_model)
```

```
WARNING:absl:Found untraced functions such as conv2d_layer_call_fn, conv2d_layer_call_ar
```

```
WARNING:absl:Found untraced functions such as conv2d_layer_call_fn, conv2d_layer_call_ar
```

```
INFO:tensorflow:Assets written to: /tmp/tmpvjcj690u/assets
```

```
INFO:tensorflow:Assets written to: /tmp/tmpvjcj690u/assets
```

```
7588272
```

```
interpreter_predictor = tf.lite.Interpreter(model_path=str(tflite_model_quant_file))
```

```
interpreter_predictor.allocate_tensors()
```

```
!unzip pneumothorax_classifier.zip
```

```
Archive: pneumothorax_classifier.zip
```

```
creating: pneumothorax_classifier/assets/
```

```
inflating: pneumothorax_classifier/saved_model.pb
```

```
creating: pneumothorax_classifier/variables/
```

```
inflating: pneumothorax_classifier/variables/variables.data-00000-of-00001
```

```
inflating: pneumothorax_classifier/variables/variables.index
```

```
pneumothorax_classifier = tf.keras.models.load_model("pneumothorax_classifier")
```

```
#Classifier
```

```
converter_class = tf.lite.TFLiteConverter.from_keras_model(pneumothorax_classifier)
```

```
tflite_models_dir = pathlib.Path("/content/pneumothorax/")
```

```
tflite_models_dir.mkdir(exist_ok=True, parents=True)
```

```
converter_class.optimizations = [tf.lite.Optimize.DEFAULT]
```

```
tflite_classifier_model = converter_class.convert()
```

```
tflite_model_quant_file = tflite_models_dir/"classifier.tflite"
```

```
tflite_model_quant_file.write_bytes(tflite_classifier_model)
```

```
WARNING:absl:Found untraced functions such as conv2d_36_layer_call_fn, conv2d_36_layer_c
```

```
WARNING:absl:Found untraced functions such as conv2d_36_layer_call_fn, conv2d_36_layer_c
```

```
INFO:tensorflow:Assets written to: /tmp/tmp1_l3891a/assets
```

```
INFO:tensorflow:Assets written to: /tmp/tmp1_l3891a/assets
```

```
6796272
```

```
interpreter_classifier = tf.lite.Interpreter(model_path=str(tflite_model_quant_file))
```

```
interpreter_classifier = tf.lite.Interpreter(model_path=str(tflite_model_quant_file))
interpreter_classifier.allocate_tensors()
```

```
!unzip pneumothorax_classifier_oversampling.zip
```

```
Archive:  pneumothorax_classifier_oversampling.zip
  creating: pneumothorax_classifier_oversampling/assets/
  inflating: pneumothorax_classifier_oversampling/saved_model.pb
  creating: pneumothorax_classifier_oversampling/variables/
  inflating: pneumothorax_classifier_oversampling/variables/variables.data-00000-of-0006
  inflating: pneumothorax_classifier_oversampling/variables/variables.index
```

```
pneumothorax_classifier_oversampling = tf.keras.models.load_model("pneumothorax_classifier_ov
```

```
#Classifier_Oversampled
```

```
converter_class_oversampled = tf.lite.TFLiteConverter.from_keras_model(pneumothorax_classifie
tflite_models_dir = pathlib.Path("/content/pneumothorax/")
tflite_models_dir.mkdir(exist_ok=True, parents=True)
converter_class.optimizations = [tf.lite.Optimize.DEFAULT]
tflite_classifier_model = converter_class_oversampled.convert()
tflite_model_quant_file = tflite_models_dir/"classifier_over.tflite"
tflite_model_quant_file.write_bytes(tflite_classifier_model)
```

```
WARNING:absl:Found untraced functions such as conv2d_36_layer_call_fn, conv2d_36_layer_c
WARNING:absl:Found untraced functions such as conv2d_36_layer_call_fn, conv2d_36_layer_c
INFO:tensorflow:Assets written to: /tmp/tmpghslbavs/assets
INFO:tensorflow:Assets written to: /tmp/tmpghslbavs/assets
26403360
```

```
interpreter_classifier_oversampled = tf.lite.Interpreter(model_path=str(tflite_model_quant_fil
interpreter_classifier_oversampled.allocate_tensors())
```

```
#Sampling 200 images because of resource constraints
train_sampled = train.sample(200)
```

```
#Getting the overall metric of model non - oversampled
```

```
dice_coef_avg = 0
for index in tqdm(range(len(train_sampled))):
    image = tf.keras.preprocessing.image.load_img(train_sampled.iloc[index]['images_paths'])
    image = tf.keras.preprocessing.image.img_to_array(image, dtype='float32')
    image = tf.image.resize(image, [img_size, img_size])
    image = image / 255.0
    image = exposure.equalize_adapthist(image)      # contrast correction

    mask = tf.keras.preprocessing.image.load_img(train_sampled.iloc[index]['mask'], color_mod
    mask = tf.keras.preprocessing.image.img_to_array(mask, dtype = 'float32')
    mask = tf.image.resize(mask, [img_size, img_size])
    mask = mask / 255.0
```

```
mask = mask / 255.0
```

```
image = tf.expand_dims(image, axis = 0)
```

```
is_pneumothorax = pneumothorax_classifier.predict(image)
```

```
is_pneumothorax = is_pneumothorax[0]
```

```
if is_pneumothorax >= threshold :
```

```
    pred_masks = pneumothorax_predictor.predict(image)
```

```
    pred_masks = pred_masks[0]
```

```
else :
```

```
    pred_masks = np.zeros((256,256,1))
```

```
dice_coef_avg = dice_coef_avg + dice_coef(mask, pred_masks)
```

```
100%|██████████| 200/200 [01:28<00:00, 2.25it/s]
```

```
print("Average dice coef for non - oversampled model is:", dice_coef_avg.numpy()/200)
```

```
Average dice coef for non - oversampled model is: 0.11101271629333496
```

```
#Getting the overall metric of quantized model non - oversampled
```

```
dice_coef_avg = 0
```

```
for index in tqdm(range(len(train_sampled))):
```

```
    image = tf.keras.preprocessing.image.load_img(train_sampled.iloc[index]['images_paths'])
```

```
    image = tf.keras.preprocessing.image.img_to_array(image, dtype='float32')
```

```
    image = tf.image.resize(image, [img_size, img_size])
```

```
    image = image / 255.0
```

```
    image = exposure.equalize_adapthist(image)      # contrast correction
```

```
    mask = tf.keras.preprocessing.image.load_img(train_sampled.iloc[index]['mask'], color_mod
```

```
    mask = tf.keras.preprocessing.image.img_to_array(mask, dtype = 'float32')
```

```
    mask = tf.image.resize(mask, [img_size, img_size])
```

```
    mask = mask / 255.0
```

```
image = tf.expand_dims(image, axis = 0)
```

```
input_index = interpreter_classifier.get_input_details()[0]["index"]
```

```
output_index = interpreter_classifier.get_output_details()[0]["index"]
```

```
interpreter_classifier.set_tensor(input_index, image)
```

```
interpreter_classifier.invoke()
```

```
is_pneumothorax = interpreter_classifier.get_tensor(output_index)
```

```
is_pneumothorax = is_pneumothorax[0]
```

```
if is_pneumothorax >= threshold :
```

```
    input_index = interpreter_predictor.get_input_details()[0]["index"]
```

```
    output_index = interpreter_predictor.get_output_details()[0]["index"]
```

```
    interpreter_predictor.set_tensor(input_index, image)
```

```
    interpreter_predictor.invoke()
```

```
    pred_masks = tf.reshape(interpreter_predictor.get_tensor(output_index), [256, 256, 1])
```



```

else :
    pred_masks = np.zeros((256,256,1))

dice_coef_avg = dice_coef_avg + dice_coef(mask, pred_masks)

100%|██████████| 200/200 [44:39<00:00, 13.40s/it]

print("Average dice coef for quantized non - oversampled model is:", dice_coef_avg.numpy()/20

    Average dice coef for quantized non - oversampled model is: 0.07591805458068848

#Getting the overall metric of model oversampled
dice_coef_avg = 0
for index in tqdm(range(len(train_sampled))):
    image = tf.keras.preprocessing.image.load_img(train_sampled.iloc[index]['images_paths'])
    image = tf.keras.preprocessing.image.img_to_array(image, dtype='float32')
    image = tf.image.resize(image, [img_size, img_size])
    image = image / 255.0
    image = exposure.equalize_adapthist(image)      # contrast correction

    mask = tf.keras.preprocessing.image.load_img(train_sampled.iloc[index]['mask'], color_mod
    mask = tf.keras.preprocessing.image.img_to_array(mask, dtype = 'float32')
    mask = tf.image.resize(mask, [img_size, img_size])
    mask = mask / 255.0

    image = tf.expand_dims(image, axis = 0)

    is_pneumothorax = pneumothorax_classifier_oversampling.predict(image)
    is_pneumothorax = is_pneumothorax[0]

    if is_pneumothorax >= threshold :
        pred_masks = pneumothorax_predictor.predict(image)
        pred_masks = pred_masks[0]
    else :
        pred_masks = np.zeros((256,256,1))

    dice_coef_avg = dice_coef_avg + dice_coef(mask, pred_masks)

100%|██████████| 200/200 [01:48<00:00, 1.84it/s]

print("Average dice coef for oversampled model is:", dice_coef_avg.numpy()/200)

    Average dice coef for oversampled model is: 0.142722749710083

```

```

#Getting the overall metric of quantized model oversampled
dice_coef_avg = 0
for index in tqdm(range(len(train_sampled))):
    image = tf.keras.preprocessing.image.load_img(train_sampled.iloc[index]['images_paths'])
    image = tf.keras.preprocessing.image.img_to_array(image, dtype='float32')
    image = tf.image.resize(image, [img_size, img_size])

```

```

image = tf.image.resize(image, [img_size, img_size])
image = image / 255.0
image = exposure.equalize_adapthist(image)      # contrast correction

mask = tf.keras.preprocessing.image.load_img(train_sampled.iloc[index]['mask'], color_mod
mask = tf.keras.preprocessing.image.img_to_array(mask, dtype = 'float32')
mask = tf.image.resize(mask, [img_size, img_size])
mask = mask / 255.0

image = tf.expand_dims(image, axis = 0)

input_index = interpreter_classifier_oversampled.get_input_details()[0]["index"]
output_index = interpreter_classifier_oversampled.get_output_details()[0]["index"]
interpreter_classifier_oversampled.set_tensor(input_index, image)
interpreter_classifier_oversampled.invoke()
is_pneumothorax = interpreter_classifier_oversampled.get_tensor(output_index)
is_pneumothorax = is_pneumothorax[0]

if is_pneumothorax >= threshold :
    input_index = interpreter_predictor.get_input_details()[0]["index"]
    output_index = interpreter_predictor.get_output_details()[0]["index"]
    interpreter_predictor.set_tensor(input_index, image)
    interpreter_predictor.invoke()
    pred_masks = tf.reshape(interpreter_predictor.get_tensor(output_index), [256, 256, 1])

else :
    pred_masks = np.zeros((256,256,1))

dice_coef_avg = dice_coef_avg + dice_coef(mask, pred_masks)

100%|██████████| 200/200 [43:23<00:00, 13.02s/it]

print("Average dice coef for quantized oversampled model is:", dice_coef_avg.numpy()/200)

Average dice coef for quantized oversampled model is: 0.09985686302185058

def display(display_list):
    plt.figure(figsize=(25, 25))

    title = ['Input Image', 'True Mask', 'Predicted Mask', 'Predicted Mask(Quantized)', 'Pred

    for i in range(len(display_list)):
        plt.subplot(1, len(display_list), i+1)
        plt.title(title[i])
        plt.imshow(tf.keras.preprocessing.image.array_to_img(display_list[i]))
        plt.axis('off')
    plt.show()

train_display = train.sample(5)

```

```

for index in range(len(train_display)):
    image = tf.keras.preprocessing.image.load_img(train_display.iloc[index]['images_paths'])
    image = tf.keras.preprocessing.image.img_to_array(image, dtype='float32')
    image = tf.image.resize(image, [img_size, img_size])
    image = image / 255.0
    image = exposure.equalize_adapthist(image)      # contrast correction

    mask = tf.keras.preprocessing.image.load_img(train_display.iloc[index]['mask'], color_mod
    mask = tf.keras.preprocessing.image.img_to_array(mask, dtype = 'float32')
    mask = tf.image.resize(mask, [img_size, img_size])
    mask = mask / 255.0

    image = tf.expand_dims(image, axis = 0)

    #Non - Oversampled
    is_pneumothorax = pneumothorax_classifier(image)
    is_pneumothorax = is_pneumothorax[0]

    if is_pneumothorax >= threshold :
        pred_masks = pneumothorax_predictor.predict(image)
        pred_masks = pred_masks[0]
    else :
        pred_masks = np.zeros((256,256,1))

    #Quantized
    input_index = interpreter_classifier.get_input_details()[0]["index"]
    output_index = interpreter_classifier.get_output_details()[0]["index"]
    interpreter_classifier.set_tensor(input_index, image)
    interpreter_classifier.invoke()
    is_pneumothorax = interpreter_classifier.get_tensor(output_index)
    is_pneumothorax = is_pneumothorax[0]

    if is_pneumothorax >= threshold :
        input_index = interpreter_predictor.get_input_details()[0]["index"]
        output_index = interpreter_predictor.get_output_details()[0]["index"]
        interpreter_predictor.set_tensor(input_index, image)
        interpreter_predictor.invoke()
        pred_masks_quant = tf.reshape(interpreter_predictor.get_tensor(output_index), [256, 2
    else :
        pred_masks_quant = np.zeros((256,256,1))

    #Oversampled
    is_pneumothorax = pneumothorax_classifier_oversampling(image)
    is_pneumothorax = is_pneumothorax[0]


    if is_pneumothorax >= threshold :
        pred_masks_over = pneumothorax_predictor.predict(image)
        pred_masks_over = pred_masks_over[0]
    else :
        pred_masks_over = np.zeros((256,256,1))

```

```
#Quantized Oversampled
input_index = interpreter_classifier_oversampled.get_input_details()[0]["index"]
output_index = interpreter_classifier_oversampled.get_output_details()[0]["index"]
interpreter_classifier_oversampled.set_tensor(input_index, image)
interpreter_classifier_oversampled.invoke()
is_pneumothorax = interpreter_classifier_oversampled.get_tensor(output_index)
is_pneumothorax = is_pneumothorax[0]

if is_pneumothorax >= threshold :
    input_index = interpreter_predictor.get_input_details()[0]["index"]
    output_index = interpreter_predictor.get_output_details()[0]["index"]
    interpreter_predictor.set_tensor(input_index, image)
    interpreter_predictor.invoke()
    pred_masks_quant_over = tf.reshape(interpreter_predictor.get_tensor(output_index), [2
else :
    pred_masks_quant_over = np.zeros((256,256,1))

display([image[0], mask, pred_masks, pred_masks_quant, pred_masks_over , pred_masks_quant
```



-- Quantized models have reduced quality

